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00100EX-20161026

0483

Vario-Soft-Geschiebe

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(EN) Instructions for use

Vario-Soft attachments

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non compliance with the instructions for use.

1. Description of the product

The Vario-Soft attachments are rod attachments with friction matrices in green (reduced friction, 4 N), yellow (medium / normal friction, 6 N) and red (strong friction, 8 N) for the fabrication of removable dental restorations. Color-coding of the matrices enables fast and accurate determination of the friction. The various attachments are distinguished by their shape and size. The patrices are made of burn-out resin (polystyrene, PS); the matrices are made of polyoxymethylene (POM).

2. Indication

The Vario-Soft attachments are used for fixation of a removable restoration with friction. Depending on the required level of friction, users can choose between three color-coded matrices. For correct placement/insertion of the matrices, reference is made to the original accessories. The patrices of the attachments Vario-Soft sv (vs3 sv) and Vario-Soft 3 mini sv (vs3 mini sv) feature an integrated shear distributor. The fabrication of an individually milled shear distributor with interlock is not required when using these two attachments.

The individually milled shear distributor with interlock is always required when using the attachments Vario-Soft 3 (vs3) and Vario-Soft 3 mini (vs3 mini). The attachments are only designed for retention of the denture and not for distribution of the shear force to the fixed restoration.

Full-cover crowns, implant crowns, bar restorations, interdental gaps, free-end dentures.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. Safety hand hazard information

The lot numbers of all components should be recorded to ensure traceability and to handle complaints. Adequate protection against swallowing and aspiration must be ensured when inserting the restoration.

5. Storage and durability

The product is supplied in the non-sterile condition and dry, dust-free storage in the original packaging must be ensured.

Storage:

Please observe the storage information printed on the label.

Durability:

There are no restrictions with regard to durability.

6. Processing

6.1 Preparation

The size of the attachment needs to be included into the prosthetic planning before selecting the type of attachment. Possibility of shortening the attachments, see item 7 Technical data.

ATTENTION		
Since the Vario-Soft attachments are parallel attachments, parallel alignment of two or more patrices to one another must always be ensured to enable smooth and trouble-free insertion and removal of the denture.		

Always use the paralleling mandrel universal (REF 360 0115 1) to attach the patrices to the model using wax, resin or adhesive. Plane support on the gingiva must be ensured. If required, the plastic matrices are adapted using a tungsten carbide bur.

ATTENTION		
When using the vs3 and vs3 mini attachments, an additional individually milled shear distributor with interlock need to be attached to the model - this is not required for bar restorations.		

The primary model can be cast in the alloy of your choice.

6.2 Milling the practices

Devest the primary structure after casting, sandblast with abrasive beads and fit. After fitting and try-in in the patient's mouth, the attachment patrices and the shear distributor with interlock are finished (milled) on an individual milling base in the milling device. Milling tools with a diameter of 1 mm (e. g. REF F137 2H 10) can be used for finishing. Finishing does not result in a change to the shape of the patrices but is only a type of surface treatment to remove minor roughness. Then the patrices are polished to high gloss.

6.3 Duplicating

ATTENTION		
Polishing residues and residual milling oil are removed after polishing. Make sure to remove all residues from the cast objects.		

Preparing for duplicating is carried out based on the design of the cast structure. For fabrication of a precision-fit matrix housing, the white duplicating matrix is placed on the patrix and adapted to the gingiva if required. In the occlusal area, the duplicating matrix is fitted flush with the patrix and covered with a thin coat of wax to prevent the penetration of duplicating silicone between the patrix and duplicating matrix. Duplicating is carried out in a suitable flask using a corresponding duplicating silicone.

6.4 Investment model

The duplicating mould is filled with an investment material in the usual way. After removing the investment mold from the mould, the restoration is modelled. The wax matrix housings, which can be exactly fitted on the investment mold with the attachment matrices, allow achieving the required layer thickness of the matrix housing. This way a model with a space saving design is fabricated. Additionally, retentions can be attached later on. Sprues are attached to the model in accordance with the casting technique and investment material is filled in.

6.5 Cast metal structure

ATTENTION		
After devesting, the area of the attachment/s of the cast metal structure is sandblasted (max. 110 µm) using a sandblasting device and the investment material is removed.		

If there are any rough spots in the matrix housing, they should be thoroughly removed with a suitable too. The cast metal structure is finished and polished.

ATTENTION		
Precise fit of the individually milled or integrated shear distributor must be ensured. Correct fit and adequate function of the shear distributor are not guaranteed if the primary structure can be moved in the cast metal structure.		

6.6 Inserting the matrices

ATTENTION		
The original accessories should be used for proper insertion of the matrices.		

The matrices are placed in the housing using the corresponding insertion instrument and shortened with a scalpel if required. It is possible to choose between the various matrices in accordance with tha patient's wish.

The cast metal restoration can be completed. If the size of the matrix housing is too large and the matrices are no safely retained, they can be safely glued in place using the matrix adhesive (REF 540 0103 1). Observe the instructions for use.

6.7 Exchanging the matrices

If adequate retention of the denture is not ensured or if the patient has problems when removing his denture, the matrices can be exchanged. For this purpose, use an instrument to remove the matrices from the housing, clean the housing and insert the corresponding matrices with the insertion pin. Use a scalpel to adapt to the shape of the cast metal structure from the basal direction.

6.8 Vario-Soft 3 matrix housing

The prefabricated matrix housing is used for fabrication later on or for stress-free bonding/luting to the metal matrix housing for the Vario-Soft 3 and Vario-Soft 3 sv attachments.

ATTENTION		
A shear distributor is always required for the matrix housing.		

A yellow matrix is placed into the matrix housing and placed on the primary structure. The appendix and the housing can be customized. Then the shear distributor is modelled and the housing is cast in the alloy of your choice. The housing with the shear distributor is fitted on after casting. For the fabrication of the cast metal structure (restoration), the duplicating matrix is placed on the appendix and the model is duplicated. The restoration is waxed up and cast.

After finishing the cast metal structure, the inner surface of the housing is sandblasted and bonded in a way to eliminate stress using DTK adhesive. Observe the instructions for use. The cast metal structure (restoration) can be completed.

7. Technical data

The patrices (blue) are made of polystyrene (PS). The white, green, yellow and red matrices are made of polyoxyme-thylene (POM).

Additionally softer matrices in light green, light yellow and light green made of polypropylene (PP) are available for the Vario-Soft 3 and Vario-Soft 3 sv attachments.

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8. Symbols

CE	European conformity marking
	Caution: Federal (USA) law restricts this product to sale by or on the order of a dentist or physician.
	Reference number
	Lot number
	Manufacturer
	Attention!
	Non-sterile
	Single use
	Dry storage

9. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff.

Only original tools and components may be used for processing.

The lot numbers of all components should be recorded to ensure traceability and to handle complaints.

All products labelled are intended for single use only. However, if this product is used several times, the following risks may arise: cross contamination, malfunction, inaccurate fit, etc.

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

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