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Novo-Grip

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Instructions for use

Novo-Grip

Before using the product, please read these instructions for careful use.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.


1. Description of the product

Novo-Grip is a pair of pliers for telescopic and conical crowns and features exchangeable diamond inserts. An Allen key is used to exchange the inserts. Depending on the indication or use, inserts with different diameters are available.

Note: If excess pressure is exerted, the inserts may be bent or broken; therefore always exert adequate pressure.

2. Indication range

In the field of telescopic and conical crown techniques, Novo-Grip is used to pull primary crowns out of the secondary crowns.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Safety recommendations and hazard warnings

If excess pressure is exerted, the inserts may be bent or broken; therefore always exert adequate pressure.

5. Storage and durability

The Novo-Grip pliers and the inserts must be stored under dry and dust-free conditions.

6. Processing

6.1 Preparing Novo-Grip

Novo-Grip is supplied with two Novo-Grip inserts, two Allen screws, and an Allen key. Prior to the first use, insert the Novo-Grip inserts with the shaft side as far as possible into the bores of the Novo-Grip pliers. After that the Novo-Grip inserts slightly in order to adjust it to the same length. Place the two Allen screws onto the Allen key with the hexagon socket screw into the thread turns and tighten gently.

6.2 Use of Novo-Grip

Insert the Novo-Grip with the Novo-Grip inserts into the primary crowns and press against the inner surfaces of the crown exerting only slight pressure. Exert only slight pressure onto the handle of the pliers since the Novo-Grip inserts provide excellent grip even at minimum pressure. Exert only so much pressure to prevent the Novo-Grip inserts from sliding off. In case of excess pressure the Novo-Grip inserts may be bent or broken.

6.3 Warn Novo-Grip inserts

If the Novo-Grip inserts are worn, loosen the Allen screws and turn the Novo-Grip inserts by 90°. This way new surfaces with high grip of the diamond-coated tips come into contact at the inner surface of the crown. Roughen entirely worn Novo-Grip inserts using the Diabolo cleaner. For this purpose loosen the Allen screws, remove the Novo-Grip inserts from the pliers and clamp into the technician's handpiece. The diameter of the clamping chuck requires a diameter of 2.35 mm for this purpose. Hold the Novo-Grip inserts against the Diabolo Cleaner for 5 seconds at speed of approx. 3000 revolutions per minute and exert slight pressure. Wear of the Diabolo Cleaner is desired and absolutely normal. Then reinsert the Novo-Grip inserts into the bores of the Novo-Grip and fasten the Allen screws.

7. Accessories

Novo-Grip inserts, standard, 2 pieces	REF 310 0001 A
Novo-Grip inserts, small, 2 pieces	REF 310 0001 B
Threaded pins M3, 4 pieces	REF 310 0011 2
Diabolo-Cleaner, 1 piece	REF 340 0100 0

8. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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