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000160EX-20161110

**CE**  
0483

## Vario-Kugel-Snap Geschiebe oc

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## (EN) Instructions for use

## Vario-Stud-Snap attachments oc

**Before using the product, please read these instructions for use carefully!**

**The manufacturer will not accept any liability for damage resulting from non compliance with the instructions for use.**

1. **Description of the product**

The Vario-Stud-Snap attachments oc (vks-oc) are available in different sizes (1.7 and 2.2 mm) with snap matrices in green (reduced snap, 4 N), yellow (medium / normal friction, 6 N) and red (strong snap, 8 N) and are used for the fabrication of removable dental restorations. Color-coding of the matrices enables fast and accurate determination of the snap type. The different attachments are distinguished by their shape and size. The patrices are made of burn-out resin (polystyrene, PS) for fabrication in the alloy of your choice; the matrices are made of polyoxymethylene (POM).

2. **Indication**

The Vario-Stud-Snap attachments are used for fixation of a removable dental restoration with snap retention. Depending on the required level of snap, users can choose between three color-coded matrices. For correct placement/insertion, reference is made to the original accessories.

An individually milled shear distributor with interlock is always required when using the attachments vks-oc 1.7 30°, vks-oc 2.2 60° and the extracoronal attachment vks-oc 2.2. The attachments are only designed for retention of the denture and not for distribution of the shear force to the fixed restoration.

The patrices vks-oc/srg 1.7 or 2.2 are used for root caps and bar structures.

Full-cover crowns, implant crowns, bar restorations, interdental gaps, free-end dentures, root caps.

3. **Contraindication**

Based on the current state of knowledge there are no known contraindications.

The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. **Safety and hazard information**

The lot numbers of all components should be recorded to ensure traceability and to handle complaints. Adequate protection against swallowing and aspiration must be ensured when inserting the restoration.

5. **Storage and durability**

The product is supplied in the non-sterile condition and dry, dust-free storage in the original packaging must be ensured.

**Storage:**

Please observe the storage information printed on the label.

**Durability:**

There are no restrictions with regard to durability.

6. **Processing**

6.1 **Preparation**

The size of the attachment needs to be included into the prosthetic planning before selecting the type of attachment. Information about sizes can be found under item 7. Technical data.

<span>⚠</span> <b>ATTENTION</b>
<b>It is also essential for stud attachments to ensure parallel alignment of two or more patrices to one another to enable smooth and trouble-free insertion and removal of the denture.</b>

Always use the paralleling mandrel oc/srg REF 430 0677 0 for size 1.7 and REF 360 0113 0 for size 2.2 to attach the patrices to the model using wax resin or adhesive.

<span>⚠</span> <b>ATTENTION</b>
<b>When using the attachments vks-oc 1.7 30°, vks-oc 2.2 60° or the extracoronal attachment vks-oc 2.2, an additional individually milled shear distributor with interlock needs to be attached to the model - this is not required for bar restorations.</b>

The primary model can be cast in the alloy of your choice.

6.2 **Milling the patrices**

Devast the primary structure after casting, sandblast with abrasive beads and fit. After fitting and try-in in the patient's mouth, the attachment patrices and the shear distributors with interlock are finished (milled) and the individual milling base in the milling device. Milling tools with a diameter of 1 mm (e.g. REF F137 2H 10) can be used for finishing. Finishing does not result in a change to the shape of the patrices but is only a type of surface treatment to remove minor roughness. Then the patrices are polished to high gloss.

6.3 **Duplicating**

<span>⚠</span> <b>ATTENTION</b>
<b>Polishing residues and residual milling oil are removed after polishing. Make sure to remove all residues from the cast objects.</b>

Preparing for duplicating is carried out based on the design of the cast structure. The yellow matrix is placed onto the patrix and the blocking out disc is placed below the matrix for the fabrication of a precision-fit matrix housing. Parallel blocking out (corresponding size) of the blocking out disc towards the model is carried out. Duplicating is carried out in a suitable flask using a corresponding duplicating silicone.

6.4 **Investment model**

The duplicating mould is filled with an investment material in the usual way. After removing the investment mold from the mould, the restoration is modelled. The moulded matrix is covered with a thin coat of wax (at least 0.2 mm). Additionally, retentions can be attached later on.

Sprues are attached to the model in accordance with the casting technique and investment material is filled in.

6.5 **Cast metal structure**

<span>⚠</span> <b>ATTENTION</b>
<b>After divesting, the area of the attachment/s of the cast metal structure is sandblasted (max. 110 µm) using a sandblasting device and the investment material is removed.</b>

If there are any rough spots in the matrix housing, they should be thoroughly removed with a suitable tool. The cast metal structure is finished and polished.

<span>⚠</span> <b>ATTENTION</b>
<b>Precise fit of the individually milled shear distributor must be ensured. Correct fit and adequate function of the shear distributor are not guaranteed if the primary structure can be moved in the cast metal structure.</b>

### 6.6 Inserting the matrices

<span>⚠</span> <b>ATTENTION</b>
<b>The original accessories should be used for proper insertion of the matrices.</b>

The matrices are placed in the housing using the corresponding insertion instrument. It is possible to choose between the various matrices in accordance with the patient's wish. The cast metal restoration can be completed.

### 6.7 Exchanging the patrices

If adequate retention of the denture is not ensured or if the patient has problems when removing his denture, the matrices can be exchanged. For this purpose, use an instrument to remove the matrices from the housing, clean the housing and insert the corresponding matrices with the insertion pin.

6.8 **Root caps**

The patrix is fabricated as described above. If a restoration is fabricated without a cast substructure, the matrices can be fixed with the matrix housing directly in the acrylic resin.

The blocking out disc needs to be placed below the stud. The insertion pin is used to insert the yellow matrix in the metal patrix housing and the matrix, which is filled with liquid silicone, is pressed onto the patrix and blocked out parallelly below the blocking out disc. Then the denture can be completed.

6.9 **Repair**

If repair is required, the metal transfer patrices are inserted into the matrices to ensure reliable transfer or the position to the model. Then the denture can be repaired or relined.

7. **Technical data**

The patrices (blue) are made of polystyrene (PS). The green, yellow and red matrices are made of polyoxymethylene (POM).

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8. **Symbols**

<span>CE</span>	European conformity marking
<span>Rx only</span>	Caution: Federal (USA) law restricts this product to sale by or on the order of a dentist or physician.
<span>REF</span>	Reference number
<span>LOT</span>	Attention!
<span>⚠</span>	Lot number
<span>⚠</span>	Non-sterile
<span>⚠</span>	Manufacturer
<span>ⓧ</span>	Single use
<span>☔</span>	Dry storage

9. **Additional information**

The product may only be used by dentists, dental technicians and adequately trained dental staff. Only original tools and components may be used for processing. The lot numbers of all components should be recorded to ensure traceability and to handle complaints. All products labelled 🚫 are intended for single use only. However, if this product is used several times, the following risks may arise: cross contamination, malfunction, inaccurate fit, etc. These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

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