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## (EN) Instructions for use

## Tool set for individual screw connections M1.4/ M1.6

**Please read these instructions for use carefully before using the product.**
**The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.**

### 1. Product description

The tool set for individual screw connections is used for conditionally removable prosthetics and serves to restore the chewing function as a connecting element for primary and secondary constructions. If necessary, the screw can be removed and the removable part of the dentures cleaned.

The 10-piece tool set is used to produce individual occlusal and horizontal screw connections. The conical screw head prevents the screw connection from loosening. The countersinking of the thread portion of the screw 0.3 mm into the primary element enables the shear forces generated to be transferred from the screw thread to the more stable screw head.

The product assortment includes the following components:

- One HM-centring drill
- One Diati-Multidrill
- One Diati-Multidrill with stop
- One HM facing cutter
- One first tap
- One last tap
- One tap holder
- One screwdriver, short
- Two titanium screws

### 2. Indication

Production of individual screw connections

### 2.1 Target patient groups

The patient group consists of adults with residual dentition in the upper/lower upper/lower dental dentition in com-bination with implants or patients with implant-supported restorations in the upper/lower jaw.

### 2.2 Intended users

The product must only be used by dentists, dental technicians and adequately trained dental staff.

### 3. Contraindication

The tool kit for individual screw connections is contraindicated in the following cases:

- Screw connections are not possible in CoCr alloys as these alloys have a high hardness.
- If intolerance is suspected, the product must not be used until allergy tests that confirm the absence of an allergy have been completed.

### 4. Hazard and safety information

The batch numbers of all products used must be recorded to ensure traceability and for the purpose of handling complaints.

The storage conditions (see section 5) and the handling conditions (see section 6) must be observed to ensure that the products are processed as intended. Incorrect storage can affect the product properties and result in failure of the supply.

All parts should be protected against aspiration and swallowing during use. The aspiration and swallowing of products can lead to infection and physical injury.

All non-sterile parts must be disinfected before being inserted into the patient. For this, all assembled parts must be disassembled into their individual parts.

To avoid cross-contamination the reprocessing instructions from the manufacturer must be observed.

If the secondary construction is removed for cleaning, new retention screws must be used. Reuse of the retaining screws is contraindicated.

The titanium screw is intended for single use only. If this product is nevertheless used more than once, the fol-lowing risks may arise: Cross-contamination, malfunction, poor fit etc.

### Prevention:

Read and understand all safety instructions before use. Use personal protective equipment as required.

### Product contains:

	Analysis for guide values %								
Name:	C	Fe	O	N	H	Al	V	Y	Ti
Titanium	max. 0,08	max. 0,25	max. 0,13	max. 0,05	max. 0,012	5,50–6,50	3,50–4,50	max. 0,005	Rest

Auxiliary instruments may contain nickel. A list of the materials used can be found in the table in section 7 "Tech-nical data". The medical devices for use in combination with the product can be found in the table in section 7.1.

### 5. Storage and shelf life information

The product is delivered in a non-sterile condition and must be stored in its original packaging in a dry, dust-free area and at room temperature.

### Storage:

Please observe storage instructions indicated on the label.

### Disposal:

Products, instruments and accessories that are not subject to any special requirements with regard to collection and disposal in view of the prevention of infection, may be disposed of together with household waste. The regional, national and international regulations must be observed.

### 6. Processing/handling

#### ⚠ Caution:

As a general principal, we recommend thorough surgical and prosthetic planning and selection of the prosthetic components based on the specific situation and restoration.

All parts should be protected against aspiration and swallowing during insertion. The aspiration of products can lead to infection and physical injury.

All non-sterile parts must be disinfected before being inserted into the patient. For this, all assembled parts must be disassembled into their individual parts.

The manufacturer's instructions for use of the disinfectant agent must be observed.

### 6.1 Use

The primary element must be made of gold alloy or titanium. If harder alloys are used, tapping is only possible in rare cases. To achieve an optimal result in the form of a smooth drilling process (no jamming and therefore no tilting) and a long drill service life, a lot of milling and drilling oil (REF 55000008) should be used for each work step.

#### ⚠ Caution:

The primary element must have a minimum thickness of 3.9 mm and the secondary element a minimum thickness of 0.9 mm at the point of the screw connection.

### 6.1 HM-Centring drill

The primary and secondary elements are assembled and a small groove is created using the centring drill. The groove provides adequate fixation for the Multidrill.

### 6.3 Diati-Multidrill

Use the Multidrill to drill in the defined direction through the attached secondary element approx. 1.5 to 2 mm deep into the primary element.

### 6.4 Diati-Multidrill with stop

Remove the secondary element. Insert the Multidrill with a stop into the core hole in the primary element and drill up to the stop. After drilling, remove all the metal shavings and oil residues on the primary and secondary elements. There must be no chip residues in the core hole.

### 6.5 HM facing cutter

The primary element and secondary element are assembled. Lower the HM facing cutter through the secondary element until the guide pin of the facing cutter reaches the stop in the core hole of the primary element.

#### ⚠ Caution:

The ideal speed is 5000 to 10000 rpm. Work with a copious amount of oil in this work step. The facing cutter can be clogged up by metal shavings. To maintain the facing cutter's high cutting performance, it must be cleaned with a steam jet device.

Creation of the conical countersunk in the primary and secondary elements to hold the screw head. After lowering, lift off the secondary element and insert the titanium screw into the countersunk. The screw thread must fit through the secondary element and the screw head should protrude approx. 0.3 mm from the secondary element at the thread. When the primary and secondary elements have been screwed together, the screw head is lowered 0.3 mm deep into the primary element.

### 6.6 First and last taps

Fix the first tap in the tap holder and turn it clockwise without pressure, using a lot of oil, to cut the thread into the primary element.

Procedure: Two to three turns to the right, then one turn to the left. Continue turning the tap in this way until the first tap reaches the bottom of the core hole. Unscrew the first tap by turning it to the left and clean the core hole with compressed air.

Fix the last tap in the thread tap holder and use the same procedure used with the first tap to finish cutting the thread. Thoroughly clean the core hole to remove metal chips and oil residues after cutting the thread.

### 6.7 Titanium screw

After assembling the primary and secondary construction, the titanium screw is screwed in using the screwdriver. Adjust the screw head according to the contour of the secondary element. We recommend the titanium process-ing set REF 350000890 for the grinding and polishing of the titanium screw.

### 6.8 Auxiliary modelling element

The auxiliary modelling element is placed in the core hole of the primary construction. The auxiliary modelling element is covered with the residue-free burning out resin Pi-Ku-Plast HP 36 REF 54000220. The anatomically reduced model of the secondary construction is produced according to standard dental technology rules. Using a suitable instrument (e.g. tongs), the auxiliary modelling element can be removed from the modelling with a rotary movement.

### 7. Technical data

REF:	Product:	Quantity:	Material:	Ø:	Length:	Thread:	Length/head:	Max. shortening:
33000016	Tool set for individual screw connections M1.6	10-piece						
33000600	Tool set for individual screw connections M1.4	10-piece						
<b>Accessories:</b>								
33000660	HM-Centring drill Ø1.4 for M1.4 and M1.6	1	Carbide					
33000650	Tungsten carbide facing cutter M1.4	1	Carbide					
33001159	Tungsten carbide facing cutter M1.6	1	Carbide					
33000671	First tap M1.4	1	Carbide					
3300116V	First tap M1.6	1	Carb					

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