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000259EX-20161110

CE 0483

Vario-Kugel-Snap Geschiebe sg

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(EN) Instructions for use

Vario-Stud-Snap attachments sg

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non compliance with the instructions for use.

1. Description of the product

The Vario-Stud-Snap attachments sg (vks-sg) are available in different sizes (1.7 and 2.2 mm) with snap matrices in green (reduced snap, 4 N), yellow (medium / normal friction, 6 N) and red (strong snap, 8 N) and are used for the fabrication of removable dental restorations. Color-coding of the matrices enables fast and accurate determination of the snap type. The different attachments are distinguished by their shape and size. The patrices are made of burn-out resin (polystyrene, PS) for fabrication in the alloy of your choice; the matrices are made of polyoxymethylene (POM).

2. Indication

The Vario-Stud-Snap attachments are used for fixation of a removable dental restoration with snap retention. Depending on the required level of snap, users can choose between three color-coded matrices. For correct placement/insertion, reference is made to the original accessories.

An individually milled shear distributor with interlock is always required when using the attachments vks-sg 1.7 and 2.2. This does not apply to the vks-sg 1.7 sv attachment that already features an integrated shear distributor. The attachments are only designed for retention of the denture and not for distribution of the shear force to the fixed restoration. Full-cover crowns, implant crowns, bar restorations, interdental gaps, free-end dentures.

3. Contraindication

Based on the current state of knowledge there are no known contraindications. The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. Safety and hazard information

The lot numbers of all components should be recorded to ensure traceability and to handle complaints. Adequate protection against swallowing and aspiration must be ensured when inserting the restoration.

5. Storage and durability

The product is supplied in the non-sterile condition and dry, dust-free storage in the original packaging must be ensured. **Storage:** Please observe the storage information printed on the label.

Durability:

There are no restrictions with regard to durability.

6. Processing

6.1 Preparation

The size of the attachment needs to be included into the prosthetic planning before selecting the type of attachment. Information about sizes can be found under item 7. Technical data.

	ATTENTION
	It is also essential for stud attachments to ensure parallel alignment of two or more patrices to one another to enable smooth and trouble-free insertion and removal of the denture.

Always use the paralleling mandrel oc/sjg REF 430 0677 0 for size 1.7 and REF 360 0113 0 for size 2.2 to attach the patrices to the model using wax resin or adhesive.

	ATTENTION
	When using the attachments vks-sg 1.7 and 2.2, an additional individually milled shear distributor with interlock needs to be attached to the model - this is not required for bar restorations.

The primary model can be cast in the alloy of your choice.

6.2 Milling the patrices

Devest the primary structure after casting, sandblast with abrasive beads and fit. After fitting and try-in in the patient's mouth, the attachment patrices and the shear distributors with interlock are finished (milled) on an individual milling base in the milling device. Milling tools with a diameter of 1 mm (e.g. REF F137 2H 10) can be used for finishing. Finishing does not result in a change to the shape of the patrices but is only a type of surface treatment to remove minor roughness. Then the patrices are polished to high gloss.

6.3 Duplicating

	ATTENTION
	Polishing residues and residual milling oil are removed after polishing. Make sure to remove all residues from the cast objects.

Preparing for duplicating is carried out based on the design of the cast structure. The yellow matrix is placed onto the patrix for the fabrication of a precision-fit matrix housing. Below the matrix, the space towards the gingiva is blocked out. Duplicating is carried out in a suitable flask using a corresponding duplicating silicone.

6.4 Investment model

The duplicating mould is filled with an investment material in the usual way. After removing the investment model from the mould, the restoration is modelled. The moulded matrix is covered with a thin coat of wax (at least 0.2 mm). Additionally, retentions can be attached later on.

Sprues are attached to the model in accordance with the casting technique and investment material is filled in.

6.5 Cast metal structure

	ATTENTION
	After devesting, the area of the attachment/s of the cast metal structure is sandblasted (max. 110 µm) using a sandblasting device and the investment material is removed.

If there are any rough spots in the matrix housing, they should be thoroughly removed with a suitable tool. The cast metal structure is finished and polished.

	ATTENTION
	Precise fit of the individually milled shear distributor must be ensured. Correct fit and adequate function of the shear distributor are not guaranteed if the primary structure can be moved in the cast metal structure.

6.6 Inserting the matrices

	ATTENTION
	The original accessories should be used for proper insertion of the matrices.

The matrices are placed in the housing using the corresponding insertion instrument. It is possible to choose between the various matrices in accordance with the patient's wish. The cast metal restoration can be completed.

6.7 Exchanging the matrices

If adequate retention of the denture is not ensured or if the patient has problems when removing his denture, the matrices can be exchanged. For this purpose, use an instrument to remove the matrices from the housing, clean the housing and insert the corresponding matrices with the insertion pin.

6.8 Repair

If repair is required, the metal transfer patrices are inserted into the matrices to ensure reliable transfer or the position to the model. Then the denture can be repaired or relined.

7. Technical data

The patrices (blue) are made of polystyrene (PS). The green, yellow and red matrices are made of polyoxymethylene (POM).

	vks-oc/sjg uni 1.7	vks-oc/sjg 2.2	vks-sg 1.7	vks-sg/sv 1.7	vks-sg 2.2
Height	2.2 mm	3.2 mm	4.1 mm	4.5 / 5.5 mm	5.4 mm
Width	-	-	3.0 mm	3.5 mm	3.8 mm
Depth	-	-	2.2 mm	4.3 mm	3.5 mm
Stud diameter	1.7 mm	2.2 mm	1.7 mm	1.7 mm	2.2 mm

7.2 Matrices

	vks-sg 1.7	vks-sg 2.2
Height	3.1 mm	4.1 mm
Width	3.2 mm	4.2 mm
Depth	2.3 mm	2.85 mm

7.3 Matrix housings

	vks-sg 1.7	vks-sg 2.2
Height	3.5 mm	4.5 mm
Width	4.1 mm	5.1 mm
Depth	10.5 mm	11 mm


8. Symbols

 European conformity marking

Rx only	Caution: Federal (USA) law restricts this product to sale by or on the order of a dentist or physician.
	Reference number
	Lot number
	Manufacturer
	Attention!
	Non-sterile
	Single use
	Dry storage

9. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff. Only original tools and components may be used for processing.

The lot numbers of all components should be recorded to ensure traceability and to handle complaints. All products labelled  are intended for single use only. However, if this product is used several times, the following risks may arise: cross contamination, malfunction, inaccurate fit, etc. These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

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