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6. Processing and Use 6.1 Preparation

 **CAUTION**
We recommend thorough surgical and prosthetic planning as a matter of principle and selection of the prosthetic components based on the situation and restoration.

All parts are to be secured against aspiration and swallowing during insertion. Aspiration of products can lead to infections and physical injuries.

The retention screws used in the laboratory/milling center must not be used clinically, as they are subject to loading in the laboratory, meaning that their long-term and secure positioning cannot be guaranteed.

All non-sterile parts must be prepared prior to use. All mounted parts are also to be dismantled into their individual parts.

Disinfection is then carried out using the thermo-disinfector with mildly alkaline detergents with a pH value of approx. 10,3. Sterilization of the parts can be carried out using steam sterilization via the vacuum process. To do so, a triple fractionated pre-vacuum must be generated with a sterilisation time of 4 minutes and a temperature of 134°C ± 1°C.

The manufacturer's instructions on the use of the disinfectant must be observed.

6.2 Use of the exchangeable stud vks-oc/g
Modeling primary crowns with thread sleeve
The primary model is fabricated according to standard procedures in dental techniques. In order to fix the thread sleeve to the model, a flat surface according to the largest diameter of the sleeve is required. The stud-head screw is turned in the thread sleeve and fixed to the primary crowns in the previously determined direction of insertion using the corresponding paralleling mandrel. Make sure that to remove residual wax below the stud head and the thread sleeves after completion of the model.

 CAUTION	
Long-term proper function of Vario-Stud-Snap exchangeable stud is only ensured if a lingual shear distributor with mesial interlock (inlet angle: 0°) is prepared at the abutment crown.	

Inserting and casting with thread sleeve
Prior to investing, the stud-head screw must be removed from the thread sleeve (anticlockwise) (using the short hexagon screwdriver resp. the screwdriver for the stud-head screw). To protect the thread during the casting process, a coat of colloidal graphite REF 54000706 is applied to the threaded part of the fixation screw, which is screwed in the thread sleeve subsequently. Then the model including the fixation screw is invested and cast.

 CAUTION	
When using the thread sleeve HL (melting range: 1320° C to 1460° C), the temperature of the alloy to be cast on must not exceed 1270° C.	

CAUTION
When using the platinum-iridium thread sleeve (melting range: 1820° C to 1850° C), all commercially available dental alloys except of titanium alloys can be used.

Devesting and casting
The primary structure is devested after casting and residual investment material is carefully removed with abrasive beads. The fixation screw is removed, the framework is finished and cleaned.

The suitable screwdriver is used to secure the titanium stud-head screw into the clean thread and the entire metal structure is polished to high gloss using titanium polishing paste.

 CAUTION	
Make sure that extended sandblasting does not cause any dimensional changes in the treated area.	

Modeling primary crowns with auxiliary modeling element
The primary model is fabricated according to standard methods used in dental techniques. The auxiliary modeling element is placed to the framework model using the paralleling mandrel vks-oc/g and coated with wax in the required position. Thanks to the integrated stud of the auxiliary element, the definitive position of the stud head can already be checked during the modelling process. The auxiliary modeling element is carefully removed (with a single turn) from the model after the casting process.

The auxiliary modeling element serves as a spacer for the titanium thread sleeve which is glued in the completed framework later on. (See item gluing in the thread sleeve)

Gluing in the titanium thread sleeve
The glue points are sandblasted and residual binding material is removed. The stud-head screw is screwed in the titanium thread sleeve glued in the correct position in the primary structure using DTK-adhesive REF 5400106 and the paralleling mandrel. Processing of DTK-adhesive is described in the corresponding instructions for use.

 CAUTION	
The stud-head screw in the thread sleeve using a suitable adhesive before the primary structure is inserted into the patient's mouth. Make sure not to use DTK-adhesive for fixation of the screw.	

Use of the bar patrx vks-vg titanium
Placing and adjusting

 CAUTION	
Parallel alignment of two or more bar patrices to one another must always be ensured to enable smooth and trouble-free insertion and removal of the denture.	

The titanium bar patrx vks-vg 1/7/2 exchangeable stud is shortened according to the size of the gap and the path of insertion and fitted between the crowns or implant abutments using the paralleling mandrel REF 43002860. Use the rotary instruments for titanium processing REF 35000890 to cut the bar patrx to the suitable length.

 CAUTION	
It must be ensured not to reduce the thickness of the bar patrx by milling or grinding.	

Connecting
Laser welding is used to connect the bar patrx with the crowns or the implant abutments. The operating manual of the laser welding device in use must be adhered to.

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