## Doppel-T-Klebeverbindung DTK

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## Double-T Adhesive Connector DTK

### Double-T Adhesive Connector DTK mini

Before using this product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

#### 1. Description of the product

Double-T adhesive connectors are auxiliary wax components consisting of a patrix and a matrix for simple, fast and stress-free connections between telescopic or conical crowns and CoCr frameworks. The various types ensure a wide indication range.

#### 2. Indication

The double-T adhesive connector is used for telescopic and conical crowns to ensure stress-free gluing to the CoCr framework. The double-T adhesive connector can be attached to the telescopic resp. conical crowns first as well as to the CoCr framework.

Patrices with different inclination angles – depending on the inclination of the jaw ridge – are available for the double-T adhesive connector or the double-T adhesive connector mini.

Double-T adhesive connector 90° and 120°

Double-T adhesive connector mini 90°, 105° and 120°.

#### Please note:

For stability reasons size B must only be used for saddles in partial dentures or in areas that are exposed to similar minor loads. The same applies to double-T adhesive connector mini front.

The double-T adhesive connector must only be used in conjunction with a CoCr plate or a sublingual bar.

#### 3. Contraindication

Based on the current state of knowledge there are no known contraindications.

#### 4. Storage and durability

The double-T connectors must be stored under dry conditions and not exposed to direct sunlight.

#### 5. Processing (example of DTK-mini)

##### 5.1 Double-T adhesive connector to telescopic resp. conical crowns first

The secondary crowns are modelled in the usual way after completing the primary crowns. Select inclination angle of 90°, 105° or 120° depending on the inclination of the jaw ridge. The patrices are waxed to the secondary crowns with the paralleling mandrel. It is not required, however, to achieve parallelism with the primary telescopic crowns.

Once the secondary crowns have been cast and finished, the patrices must only be sandblasted. Rough spots in the casting must be removed. Prepare the model with the secondary crowns for duplicating. Block out with wax from the patrix to the basal direction. Place the red duplicating matrix onto the patrix and duplicate.

In the duplicating mould the red duplicating matrixes must be exchanged for the blue matrix and placed in the correct position. Then pour out the duplicating mould with investment material. When destressing the investment model the blue matrixes have already been firmly integrated on the model. The CoCr modellation is carried out in the usual way and the blue matrixes are integrated.

After casting, the matrixes are only sandblasted and rough spots are removed. The CoCr framework is finished and polished in the usual way.

##### 5.1.1 Glueing

Sandblast the double-T adhesive connector with aluminium oxide 110 µm at a maximum pressure of 4 bar and then clean thoroughly. Cleaned surfaces must not be touched prior to glueing.

Check the fit of the secondary crowns and the CoCr structure. Fill the matrix of the double-T adhesive connector with DTK adhesive REF 540 0010 6 and place it on the patrix. The correct position of the CoCr structure on the model must be checked. After hardening, excess adhesive is removed with tungsten carbide burs or rubber polishers.

##### 5.2 Double-T adhesive connector at the CoCr framework first

After completing the primary crowns, duplicate the model and complete the CoCr modellation. The patrices of the double-T adhesive connector are waxed with the paralleling mandrel onto the CoCr modellation in the direction of insertion of the primary crowns. A gap is prepared for the secondary element between the primary crown and the DTK patrix. Cast and finish the CoCr framework.

The secondary crowns can be modelled on the master model parallel to the fabrication of the CoCr framework. Place the CoCr framework on the master model, place the blue matrixes onto the patrices and connect the matrixes with the secondary crowns. Remove the secondary crowns, invest, cast, finish and glue to the CoCr framework as described in section 5.1.1.

#### 6. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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