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(EN) Instructions for use

Screw connection set partly assembled

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Product description

The partly assembled screw connection set is used for conditionally removable prosthetics. As a connecting element between the primary and secondary construction, it serves to restore the masticatory function. If necessary, the screw can be replaced, and the removable part of the denture can be cleaned.

The screw connection set is a partly assembled screw connection system for the production of occlusal and horizontal screw connections in dental constructions. The thread sleeves are anchored in the construction by casting on a suitable alloy.

The product range includes the following ready-made components:

- one cuff made of a semiconductor
- one thread sleeve made of a semiconductor
- one titanium screw
- one fixation screw
- one screwdriver

2. Indication

Partially assembled connecting elements for screwing together separable bridges

– High-gold or gold reduced alloys

2.1 Patient target group

The patient group consists of adults with residual teeth in the upper and lower jaw, residual teeth in combination with implants, or implant-supported restorations in the upper and lower jaw.

2.2 Intended users

The product may be used only by dentists, dental technicians, and appropriately trained specialist personnel.

3. Contraindication

The partly assembled screw connection set is contraindicated in the following cases

- when casting with non-ferrous metal alloys
- If the secondary construction is removed for cleaning, new retaining screws must be used.
- The re-use of the retaining screws is contraindicated.
- If intolerance is suspected, this product may only be used after prior allergological clarification and proof of the absence of an allergy.

4. Hazard and safety information

The lot numbers of all products used must be documented to ensure traceability and for claims. The storage (see Point 5) and handling conditions (see Point 6) must be observed to ensure that the products are processed as intended. Improper storage can affect product properties and lead to failure of the restoration. During use, all parts must be protected against aspiration and ingestion. Aspiration and ingestion of products can lead to infections and physical injuries. All non-sterile parts must be disinfected before insertion in the patient. For this, all assembled parts must be disassembled into their individual parts. To prevent cross-contamination, the reprocessing instructions of the steriliser manufacturer must be followed. The products of the screw connection set (partly assembled) are intended for single use. If this product is used several times, the following risks may arise: cross-contamination, malfunction, and inaccuracy of fit.

Prevention:
Read and understand all safety instructions before use. Use prescribed personal protective equipment.

Product contains:

Designation:	Au:	Ir:	Pd:	Pt:
Semiconductor alloys	59,5–60,5	0,1–1,2	14,5–19,5	19,5–25,4

 	Directional analysis %									
Designation:	C	Fe	O	N	H	Al	V	Y	Ti	
Titanium	Max. 0.08	Max. 0.25	Max. 0.13	Max. 0.05	Max. 0.012	5,50–6,50	3,50–4,50	Max. 0,005	Resi-dual	

Auxiliary instruments may contain nickel. The materials used are listed in the table under Point 7 Technical data.

5. Storage and shelf life information

The product is delivered in a non-sterile condition and must be stored in its original packaging dry, dust-free, and at ambient temperature.

Storage:

Please observe the storage instructions on the label.

Disposal:

Products, instruments, accessories for which the collection and disposal are not subject to any special requirements from an infection prevention point of view can be disposed of together with household waste if necessary. The regional, national, and international regulations must be observed.

6. Processing/handling

⚠	CAUTION
We always recommend thorough surgical and prosthetic planning as well as a selection of the prosthetic parts based on the situation and the restoration.	
All non-sterile parts must be disinfected before insertion in the patient. For this, all assembled parts must be disassembled into their individual parts.	
The manufacturer's instructions for the use of the disinfectant must be observed.	

6.1 Application

⚠	CAUTION
The overall length (thread sleeve, cuff, and titanium screw) is 5.6 mm. The titanium screw with the cuff can be shortened by 1.4 mm. The minimum length of the entire construction is thus 4.2 mm. A further shortening is not possible.	

6.2 Processing the thread sleeve

After modelling the primary part, an individual stud is modelled according to the insertion direction of the abutment. For stability reasons, the cast-on thread sleeve is encased in a residue-free burn-out plastic Pi-Ku-Plast HP 36 (REF 54000220/ 54000219/ 54000217/ 54000218/ 54000216) and integrated into the stud. The remaining modelling can be done with wax.

The thread sleeve (length 3.5 mm) must never be shortened – neither from the end nor from the thread. In the event of shortening, error-free operation is no longer guaranteed.

6.3 Using the fixation screw

Use the fixation screw to hold the thread sleeve in the correct position in the muffs Apply a thin layer of colloidal graphite REF 54000706 to these before screwing them into the thread sleeve. This prevents the fixation screw from burning on and allows easy removal after casting.

⚠	CAUTION
Only the fixation screw can protect the thread from the penetration of the melt. The titanium screw is not suitable as a protective element.	

6.4 Processing of the cuff

After milling, polishing, and cleaning the primary crown, the cast-on cuff is screwed onto the cast-on thread sleeve with the titanium screw. The cuff and the primary construction are coated with a residue-free burning plastic Pi-Ku-Plast HP 36 (REF 54000220/54000219/ 54000217/ 54000218/ 54000216) and carefully removed. The remaining modelling is carried out according to the usual rules in dental technology. When processing Pi-Ku-Plast HP 36 (REF 54000220/54000219/ 54000217/ 54000218/ 54000216), follow the procedure given in the instructions for use.

6.5 Processing and shortening

Cuff and titanium screw

If required, the cuff, including titanium screw, can be shortened by 1.4 mm.

The circumferential edge on the overflow ring marks the limit up to which the screw head and the overflow ring can be shortened. Both must not be shortened further. Otherwise, the hexagon socket in the screw is no longer sufficient to tighten the screw with a screwdriver. The remaining hexagon socket is still 0.6 mm after the maximum shortening.

Thread sleeve

⚠	CAUTION
The thread sleeve must never be shortened – neither from the end nor from the thread. If the length is reduced, error-free operation is no longer guaranteed.	

6.6 Embedding and casting

The cast-on overflow ring and the cast-on threaded sleeve are made of a high-gold alloy with a melting interval of 1320–1460°C. For this reason, the casting temperature of the alloy to be cast must not exceed 1350°C.

⚠	CAUTION
The melting interval of the cast-on and solderable semiconductor alloy (Au, Pd, and Pt) is 1320–1460°C. The casting temperature of the alloy to be cast on must not exceed 1350°C. Optimum results can be achieved only if this temperature is observed.	

6.7 Divesting, milling, and polishing

After casting, divest and clean as usual.

⚠	CAUTION
The fixation screw may only be removed after the cast objects have been sandblasted. Use the fixation screw only once because the function cannot be guaranteed for a second application.	

Polishing

Polishing is carried out according to the alloy.

6.8 Cleaning

The ultrasonic unit (with a bath temperature of max. 40°C) or mechanical cleaning instruments are suitable for cleaning. Cleaning solutions should be added only according to the manufacturer's instructions.

7. Technical data

REF:	Product:	Number:	Material:	Diameter:	Length:	Thread:	Length/Head:	Max. shortening:
43007351	Prefabricated screwing set	5-part						

Zubehör:									
33000700	Titanium screw M 1.4 x 2.5 mm	1	Titanium	2.1 mm	4.5 mm	M1.4 x 0.3 mm	2.5 mm	1.4 mm	
33000710	Titanium screw M1.4 x 2.5 mm	10	Titanium	2.1 mm	4.5 mm	M1.4			

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