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Silano-Pen

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(EN) Instructions for use

Silano-Pen

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

The Silano-Pen bonding system is used to create a chemical-microchemical bond between acrylate-based, light-curing of self-curing resins (polymers) and dental alloys or ceramics (incl. aluminium and zirconium oxide).

The Silano-Pen is used for surface pretreatment during the fabrication, repair or processing of crowns, Maryland bridges, brackets, CoCr frameworks, attachments, inlays, onlays, veneers and others.

The Silano-Pen bonding system includes the following products:

- Silano-Pen
- Gas cartridge
- K-Primer bonding agent
- Disposable brush
- Plastic cup
- Cleaning brushes

2. Indication

The Silano-Pen bonding system is suitable for the following types of dental work:

- Acrylic veneering of metal frameworks
- Repair of acrylic veneers
- Repair of ceramic restorations with acrylic veneering materials (composites)
- Repair of denture resins
- Conditioning of surfaces of adhesion joints
- Bonding of acrylic teeth with denture resins
- Thermal fine disinfecting and fine cleaning of surfaces

2.1 Features and properties

Highly stable and durable bonding between metal/ceramic and acrylic is based on the combination of a special gas mixture with a bonding liquid. The special gas mixture is processed with the firing device Silano-Pen. Short firing with the Silano-Pen results in fine cleaning and, simultaneously, silicate formation and activation of the surface to be processed. Subsequent application of K-Primer (please observe the separatly enclosed instructions for use) optimizes the bond between the resulting silicate layer and the acrylic structure later on.

2.2 Target patient group

no limitations

2.3 Intended user

Dentists, dental technicians and adequately trained dental staff

3. Contraindication

The Silano-Pen may not be used intraorally. The product should not be used for patients with a suspected allergy to one or more of the components of these materials until allergy tests have been completed to confirm that the patient is not allergic to the materials.

4. Safety recommendations and hazard warning

Signal word:

Danger

Hazard information:

H222-H229 Extremely flammable aerosol. Pressurized container: may burst if heated.

Safety recommendations:

Prevention:

P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P211	Do not spray on an open flame or other ignition source.
P251	Do not pierce or burn, even after use.

5. Storage and durability

Storage:

P410+P412 Protect from sunlight. Do not expose to temperatures exceeding 50 °C. Pressurized container. The shelf life printed on the packaging is guaranteed if proper storage is ensured. The materials must not be used any longer after the shelf life has expired (see label). After filling the Silano-Pen with gas and after each use, the Silano-Pen needs to be stored in a vertical position with the nozzle facing upward. This prevents clogging of the nozzle!

6. Processing

6.1 Acrylic veneering of metal frameworks with light-curing, acrylate-based resins

6.1.1 Sandblasting
The surface to be coated is sandblasted with Al₂O₃ (disposable abrasive material: grain size 110 to 150 µm) at a pressure of 3 to 4 bar and then cleaned with water - and oil-free compressed air.

6.1.2 Silano-Pen

The metal framework es evenly heated with the flame for 5 sec per cm² (see 9.1 Silano-Pen).

6.1.3 Applying the bonding agent

Once the fired surface has cooled down (below 50 °C), K-Primer is applied in accordance with the manufacturer's instructions and allowed to dry at air for approx. 30 seconds.

6.1.4 Application of opaque and dentine materials

Apply and cure opaque and dentine materials in accordance with the manufacturer's instructions.

6.2 Repair of acrylic and ceramic restorations (also denture resins)

6.2.1 Sandblasting

In order to prepare the object for new strcturing, the repair surface is sandblasted with Al₂O₃ (disposable abrasive material: grain size 110 to 150 µm) at a pressure of 3 to 4 bar and cleaned with water- and oil-free compressed air.

6.2.2 Silano-Pen

Repair surface is evenly heated with the flame.

Processing times:	Metal framework:	approx. 5 seconds/cm ²
	Ceramic:	approx. 5 seconds/cm ²
	Resin:	approx. 2 seconds/cm ²

6.2.3 Applying the bonding agent

Once the surface heated with the flame has cooled down (below 50 °C), apply K-Primer in accordance with the instructions for use and allow to dry at air for approx. 30 seconds.

6.2.4 Application of opaque, dentine or denture resin (composite)

As required, opaque, dentine or denture resin (composite) are applied and cured in accordance with the manufacturer's instructions.

6.3 Conditioning of zirconium dioxide frameworks / ceramic frameworks

Sandblast the framework, crown or bridge with Al₂O₃ (disposable abrasive material; particle size 110 µm) at a max. pressure of 2 bar. Each veneer surface is evenly heated for 5 sec with the Silano-Pen. Apply K-Primer after the heated surface has cooled down and allow to dry at air for 30 seconds. Process and cure dentine in accordance with the manufacturer's instructions.


6.4 Further possible applications

Application	Processing time Silano-Pen	Bonding agent
Integration of acrylic teeth in denture resin	Acrylic tooth approx. 2 sec/cm ² Ceramic tooth approx. 5 sec/cm ²	yes
Conditioning of DTK (adhesive) joints	patrix and matrix 5 seconds each	yes
Fine cleaning and fine disinfecting of surfaces and instruments	Acrylic 2 x 2 sec/cm ² (with break) Ceramic/metal approx. 5 sec/cm ²	no
		no

7. Operating the Silano-Pen

When delivered, the Silano-Pen is ready for operation.

Set the sliding switch ☉ of the Silano-Pen to „ON“ until you can hear that the gas mixture escapes. By pressing the spherical switch ☉ at the opposite end of the flame outlet nozzle the gas mixture is ignited.

A „pale-blue to violet“ flame can be seen. The size of the flame is adjusted with the gas flow controller ☉ (+/-). The supply of gas is interrupted by sliding back the sliding switch ☉. Operating temperature of the Silano-Pen: 15 to 25° C (also applies to filling).

8. Cleaning, filling and maintenance of the Silano-Pen

8.1 Cleaning

The flame outlet nozzle ☉ is cleaned with the enclosed small cleaning brush **each time before the Silano-Pen is filled**. For this purpose clean the nozzle with the brush for 30 seconds and push the bristles into the nozzle holes. Clean the nozzle area thoroughly with oil-free compressed air (2 to 4 bar) and remove loose particles in the nozzle openings.

8.2 Refilling with the enclosed cartridge

When refilling the Silano-Pen, the refill cartridge and the Pen must have room temperature (15 to 25 °C). Prior to filling with the enclosed refill cartridge, the flame outlet nozzle must be cleaned as described under item 8.1. Just before filling, the refill valve ☉ is shortly pressed in several times using the shaft end of the brush to ensure that the Pen is completely filled with gas. Shake the cartridge well and plug it into the refill valve of the Silano-Pen (upside down). After the cartridge has been pressed into the filling inlet, the gas mixture flows into the Silano-Pen. The filling time is 10 seconds and should not be exceeded.

8.3 Refilling with additional (max. number: 4) cartridges

Prior to using a new, closed cartridge, the sliding switch of the emptied Pen must be moved to „ON“ for at least 30 minutes to deaerate the gas tank completely. For this purpose the Silano-Pen is placed under an exhaust system or at an adequately ventilated place. Then the Pen is filled as described under item 8.2.

8.4 Checking the flame size

The length of the inner, clearly visible flame cone ☉ should be approx. 1.5 cm at room temperature and maximum position of the gas flow controller. In case of flame lengths of more than 2 cm perfect firing results can no longer be ensured. In such cases the Pen must be returned for maintenance.

Please note:

The cartridge enclosed in the set allows to prepare approx. 2000 veneer units. The period of usage of the Silano-Pen ends after a total number of 5 gas cartridges (maximum) have been filled into the Pen. According to our current knowledge, the functional reliability of the Pen is no longer ensured if the indicated period of usage is exceeded. Therefore the Pen must not be used after this period has expired.

9. Handling and processing

9.1 Silano-Pen

The surface is heated with the „pale-blue“ reactive flame zone ☉ of the pen (fig. 2). The flame is continuously moved forward and backward across the activation surface so that uniform treatment of the entire surface is achieved. The preparation of the surface to be activated and the relevant processing times are described under item 6. K-Primer is applied immediately after the heated/fired areas have cooled down below 50 °C.

Please note:

Insufficient distance or excessive distance between the flame and the surface to be treated will impair the activation result. Contamination (fingers) of the treated surfaces will affect bonding of the subsequent layers. The entire activation process, starting with sandblasting (see item 6. Processing), must be repeated in case of contamination.

9.2 K-Primer

Please observe the enclosed instructions for use when using K-Primer.

10. Symbols

	Symbol for „European Conformity“
Rx only	Caution: Federal (USA) law restricts this product to sale by