# Fräsgerät BF 2

- |    |                          |        |    |
|----|--------------------------|--------|----|
| Ⓓ  | Gebrauchsanweisung       | Seite  | 3  |
| ⒼⒷ | Instructions for use     | page   | 7  |
| Ⓕ  | Mode d'emploi            | page   | 11 |
| Ⓘ  | Istruzioni per l'uso     | pagina | 15 |
| Ⓔ  | Folleto de instrucciones | paginà | 19 |

# Fräsgerät BF 2

D

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# Milling unit BF 2

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

## 1. Intended use

The milling unit BF 2 is intended to be used in the dental laboratory. It features a three-dimensional milling arm which – as a surveyor – can be used for milling and drilling. This unit is suitable for processing materials that are common in dental laboratories, such as waxes, resins and metals.

## 2. General information

2.1 Check if voltage/power supply data are in accordance with the data on the typeplate (5.21).

2.2 The milling unit BF 2 is **not intended to be used in an environment exposed to explosion hazard**.

2.3 When using this unit, the relevant provisions of the trade association must be adhered to (working clothes and safety goggles must always be worn).

2.4 The milling unit must not be cleaned with compressed air.

2.5 To maintain the precision and extend the service life of the chuck, a tool resp. the enclosed pin must always be clamped in (even if the unit is idle).

### Please note:

- when using rotary instruments, the limits provided by the manufacturer of the instruments must be adhered to
- repair or any intervening measures may exclusively be carried out by bredent
- bredent will not accept any liability if the milling unit BF1 had not been used in accordance with these operating instructions

## 3. Installation of the unit

3.1 The unit must be placed on a level surface at working height. Safe standing must be ensured!

3.2 Provide sufficient space to move the milling arm freely.

## 4. Assembly of the unit


The milling unit needs to be assembled before it is used for the first time. For this purpose an Allen key has been enclosed to the unit.

4.1 Screw the micrometer screw (5.4) with the Allen key (4.4) and the




Allen screw to the intended position using the fork wrench.

4.2 Slide the milling spindle (4.5) into the respective locating device and tighten it with the locking screw (5.6).

4.3 The cable of the milling spindle is fed through the cable duct (5.8) and connected to the 3-pole socket (5.24) at the rear of the control unit and screwed on.

## 5. Milling unit

- (5.1) Locking screw to lock the horizontal movement
- (5.2) Locking lever for the rotation movement
- (5.3) Lever with display for drillings
- (5.4) Micrometer screw with depth stop
- (5.5) Locking screw to lock the vertical movement of the milling spindle
- (5.5) Locking screw and locating device for the milling spindle
- (5.7) Clamping lever
- (5.8) Cable bushing
- (5.9) Locking lever to fix the milling table
- (5.10) Knob for height adjustment of the milling table
- (5.11) Locking screw with height stop
- (5.12) Hole to accept the locking screw for the model platform
- (5.13) Dirt collecting tray
- (5.14) Acceptance for model or milling base
- (5.15) Locking screw
- (5.16) Locking lever for tilting and rotation movement
- (5.17) Power switch on/off
- (5.18) Right hand rotation/left-hand rotation; the milling spindle can be operated either in the left-hand or right-hand rotation mode
- (5.19) Rotary knob for continuously adjustable speed within a range of 0-30000 min<sup>-1</sup>
- (5.20) Light switch
- (5.21) Typeplate
- (5.22) Socket for power plug
- (5.23) 3-pole socket

## 6. Putting into operation

- 6.1 Check if voltage/power supply data are in accordance with the data on the typeplate (6.21).
- 6.2 Plug power plug into the socket.
- 6.3 Switch power switch on/off (5.17) to "on". **Please note:** First turn rotary knob (5.19) to "zero" to adjust the speed → **Risk of injury!**
- 6.4 Fix the model on the model platform and lock with the locking lever (5.16) in accordance with the direction of insertion.
- 6.5 Fix the model platform with the locking screw (5.12) on the milling table as required.
- 6.6 Turn the lever (5.7) to the right side to open the milling spindle, insert the desired bur and close the spindle by turning the lever to the left side. → **Tools may only be exchanged when the motor is switched off!**
- 6.7 Open the locking lever (5.9) and adjust the milling table with the knob for height adjustment (5.10) to ensure correct working height and close the locking lever.

## Milling

- 6.8.1 Open locking screws (5.1 and 5.5) and locking lever (5.2) so that milling arm and milling spindle can be moved freely. The lever (6.3) is slid to the right side out of the lock so that the milling spindle can be moved freely.
- 6.8.2 The rotary knob (5.19) is adjusted to the recommended working speed of the bur.
- 6.8.3 When milling a bar, the model is fixed on the platform to ensure parallel alignment of the milling arm and the bar. Tighten the locking lever for the rotating movement (5.2) and loosen the locking screw for horizontal movement (5.1).
- 6.8.4 Turn the rotary knob (5.19) to adjust the speed recommended for the bur.

## Drilling

- 6.9.1 Determine the desired direction of drilling at the object and close the locking screws (5.1 and 5.5) and the locking lever (5.2).
- 6.9.2 Determine the depth of the drillhole with the micrometer screw (5.4). If only the upper screw is used, the micrometer screw is not turned any longer once the depth stop has been reached. This way precise adjustment of the depth of the drillhole is ensured. The lever (5.3) is pressed in to the left until the milling spindle can be moved downward with the lever.
- 6.9.3 The rotary knob (5.19) is adjusted to the recommended working speed of the bur.
- 6.9.4 Move the milling spindle downward with the lever (5.3) whilst exerting slight pressure.


## 7. Clamping of the motor milling spindle/exchange of tools

→ Tools may only be exchanged when the motor is switched off!

- 7.1 Turn clamping lever (5.7) to the right side up to the stop and insert the tool.
- 7.2 Turn clamping lever (5.7) to the left side up to the stop and clamp in the tool.
- 7.3 By turning the clamping lever (5.7) to the right side (= approx. half of the way), the drive of the milling spindle can be locked and the tool is still firmly clamped in. This position is intended to lock a scraping tool.

## 8. Exchange of the chuck

- 9.1 Turn clamping lever (5.7) to the right side up to the stop.
- 9.2 Open locking screw of the milling spindle (5.6).
- 9.3 Remove milling spindle from the holding device.
- 9.4 Loosen union nut (8.1) completely (it may be required to use the bur shank  $\varnothing$  2.35 mm); remove (pull off) the milling spindle from the motor.
- 9.5 Turn the chuck (8.2) jerkily anticlockwise using the chuck key no. 4115 and screw it out, if required insert counter key no. 4113 in driver (8.3) for retaining purposes.
- 9.6 Apply small quantity of grease to the chuck (8.2) and insert the shaft. Using the keys, screw the chuck clockwise up to the stop as described above and tighten slightly.
- 9.7 Previously, chuck and shaft need to be carefully cleaned to avoid radial deviation.
- 9.8 Attach the milling spindle to the motor and tighten the union nut (8.1).
- 9.9 Insert milling spindle into the holding mechanism and tighten the locking screw (5.6).

## 9. Maintenance and service

### Please note:

The milling spindle must not be cleaned with compressed air!

Depending on the degree of contamination, the chuck must be cleaned and degreased (see section 8) from time to time.

## 10. Malfunctions

If the milling spindle is overloaded or blocked, the unit switches off automatically (for safety reasons). The unit can be started again by switching "off" and "on" (5.17).

If the power decreases or the motor does not run steadily, the unit should be sent to brendt.

If the spindle does not work:

- are the power supply/voltage data in accordance with the data on the typeplate?
- has the power plug been plugged in?
- has the cable of the handpiece been connected correctly to the control unit?
- power switch (5.17) must be set to "on".
- turn rotary knob (5.19) to adjust the recommended speed.
- clamping lever (5.7) must be moved to the left side.

## 11. Technical data

Width:	250 mm
Depth:	370 mm
Height:	510 mm
Weight:	19 kg
Speed:	500 - 30,000 min <sup>-1</sup>
Voltage:	100 - 115 V / 200 - 240 V
Ampere:	2 x 2.0 A
Wattage:	80 W
hertz:	50 - 60 hz

## 12. Accessories

The following accessories are available for the milling unit BF2:

Milling base	BF1 REF 140 0089 3
Tap handwheel	BF1 REF 330 0115 4
Model support	REF 730 0017 0
Chuck, 3 mm	REF 730 0015 3
Transfer device (oder: Transfer jig)	REF 360 0116 3
Brenometer surveying system	REF 310 0000 2

## 13 Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience.

The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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