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Universal Schraubendreher-Set

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(EN) Instructions for use

Universal screwdriver Set

Before using the product, please read these instructions for use carefully!

The manufacturer does not accept any liability for damage resulting from noncompliance with these instructions for use.

1. Description of the product

The universal screwdriver set includes screwdrivers that are used for the screws of various manufacturers. The screwdrivers can be used both in the laboratory and in the oral cavity. When inserting or removing a partially removable restoration or abutments, the screwdrivers are used for tightening or loosening screws. It is also possible for the dentist to activate or deactivate screws for removable dentures for friction adjustment. The screwdrivers with the torque ratchet (up to 40 Ncm) should be used for a specific torque range of the screws.

The screwdrivers are marked with a number for easy classification:

	Type	Designation	Implant system
Screwdriver 1	long/short	Torx 6	Straumann Syn-Octa
Screwdriver 2	long/short	slotted 1,6	
Screwdriver 3	long/short	slotted 2	
Screwdriver 4	short	Allen 0,03°	Altatec
Screwdriver 5	long/short	Allen 0,05°	Sulzer Calcitec-Spline Altatec Camlog Semados 3i

Screwdriver 6	long/short	Allen 0,9	Tiolox bredent
Screwdriver 7	long/short	Allen 1,0	Ankylos
Screwdriver 8	long/short	Allen 1,2	Nobel Biocare Steri-Oss Frialit II IHT Twin Plus

Screwdriver 9	long/short	Allen 1,8	Ankylos
Screwdriver 10	long	Internal Allen 2,5	bredent
Screwdriver 11	long/short	Square 1,3	3i
Screwdriver 12	long/short	Torx 5,5	Nobel Biocare

2. Indication

Replacement of lost anatomy

The screwdrivers of the Universal Screwdriver Set can be used both in the laboratory and in the oral cavity. They are used when inserting or removing a partially removable denture or abutments to tighten or loosen screws. Screws can also be activated or deactivated by the dentist for friction adjustment with removable dentures.

2.1 Patient target group

The patient group consists of adults with residual teeth in the upper/lower jaw, residual teeth in combination with implants or patients with implant-supported restorations in the upper/lower jaw.

2.2 Intended users

The intended users are dentists, dental technicians and adequately trained dental staff.

3. Contraindication

The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

- all screwdrivers that are not listed
- repeated use of screws
- unprocessed (not sterilized) and unsecured screwdrivers
- screwdrivers with damaged screw drives

4. Hazard and safety information

The lot numbers of all products used should be recorded to ensure traceability and to handle complaints. The storage and handling instructions must be observed to ensure that the products are processed as intended. Improper storage may influence the product properties and lead to failure of the restoration. All parts need to be protected against aspiration and swallowing during the use of the product. Aspiration and swallowing of products may cause infections and physical injury. In order to avoid cross contamination, the reprocessing instructions of the manufacturers of sterilization devices must be observed. Use sterile gloves or sterile instruments to remove the products to be sterilized from the packaging immediately before the use. Slipping during treatment can cause injuries that can lead to inflammation or irritation. In patients with a tendency to rapid oral aphtha formation, even minor injuries to the oral mucosa can lead to aphtha reappearance. In patients with bruxism, the use of screwdrivers needs careful consideration. On the one hand, the prosthetic restoration may be overloaded and on the other hand, the mouth opening of the patient may be severely impaired. Incorrect handling can lead to fracture of the screwdriver and aspiration or swallowing of mill components.

CAUTION

The instrument box does not serve as a sterilization box but for storage of the instruments. All parts (screwdrivers, torque ratchet and instrument box) must be sterilized individually.

Prevention

Read and understand all safety instructions before use. Use prescribed personal protective equipment.

Product contains:

The screwdrivers of the Universal Screwdriver Set contain nickel. A list of the materials used can be found in the table under item 7 Technical data.

5. Storage and durability

The product is supplied in the non-sterile condition and dry, dust-free storage at room temperature in the original packaging must be ensured. Storage after initial use only in cleaned and sterilized condition according to the instructions of the manufacturer of the sterilization device.

Storage:

Please observe the storage information on the label.

Durability:

The Universal Screwdriver Set can be reused for a maximum of 250 reprocessing cycles if the components are undamaged and not contaminated; the user is held responsible for any use beyond the number of cycles specified or the use of damaged and/or contaminated instruments.

Disposal:

Products/instruments/accessories for the collection and disposal of which no special requirements are placed as far as infection prevention is concerned can be disposed of together with household waste. Regional/national and international regulations must be observed.

Instruments are disposed of in shatter-proof, puncture-proof and sealed containers (contamination protection).

6. Processing

CAUTION

Generally, we recommend thorough surgical planning and selection of the prosthetic components depending on the situation and type of restoration. All components need to be secured against aspiration and swallowing when inserting the restoration. Aspiration of products may cause infections and physical injury. The retention screws used in the laboratory/milling center are not suitable for the clinical use, since they are subject to stress in the laboratory; as a result long-term and reliable fit can no longer be guaranteed. All non-sterile parts must be reprocessed (treated) before use (see item 6.3 Initial treatment and reprocessing). The manufacturer’s information on the use of the disinfectant must be observed.


6.1 Use

CAUTION

The screwdrivers, torque ratchet and the instrument box must be sterilized individually prior to first use in the oral cavity.

6.2 Torque ratchet

The torque ratchet needs to be disassembled in order to be cleaned and sterilized. It consists of the following parts:

Each time the ratchet is assembled and cleaned, apply maintenance oil for contra-angles or maintenance spray to the sections marked*. The torque ratchet must be calibrated after 2 years. The ratchet head and the internal spring are parts subject to wear. To insert the ratchet wheel into the ratchet head, pull the fastening pins back with thumb and forefingers, and insert the ratchet wheel at IN with the narrow side towards the ratchet head.

6.3 Initial treatment and reprocessing

Pretreatment: For treatment (cleaning/disinfection) all parts assembled must be disassembled into their individual components immediately after the use (within a maximum of 2 h) and prior to initial use and rinsed under running water for at least 1 minute. Rinse all lumens of the instruments at least three times using a disposable syringe (minimum volume 2 ml). Remove the individual parts from the instrument box. Moving parts must be moved back and forth at least three times during pre-cleaning. Place parts in pre-cleaning bath, ultrasonic cleaning (not yet activated) and brush off completely with a soft brush (move moving parts back and forth at least three times). Rinse all lumens of the instruments at least

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